

Lid-Care* Product Information

Description:	<p>Lid-Care* is a sterile, pH-adjusted, hypo-allergenic detergent solution. It contains three gentle surfactants/solubilizer (Miranol 2 MHT, PEG 20 Sorbitan monolaurate, Propylene glycol), a borate buffer is preserved with polyaminopropyl biguanide in purified water.</p> <p>Lid-Care* contains no perfumes or dyes.</p> <p>Lid-Care* Towelettes contains 4 gentle surfactants/solubilizers (PEG-200 Glyceryl Palmate, Disodium Laureth Sulfosuccinate, Coamidopropylamine Oxide, PEG-80 Glycerol cocoate).</p>
Indications:	<p>Lid-Care* is developed as a complete system for daily eyelid and lid margin hygiene. Lid-Care* is ideal as an adjunct in lid conditions associated with blepharitis. It is also an effective eye makeup remover, especially for contact lens wearers.</p>
Contraindications:	<p>Know sensitivity to any of the components</p>
Precautions:	<p>Ready to use: do not dilute. Do not purposely instill directly into the eye.</p>
Directions for use:	<ol style="list-style-type: none">1. Completely saturate a cleansing pad with Lid-Care* solution or remove a Lid-Care* Towelette from the foil pouch.2. In front of a well-illuminated mirror, fold the pad, open the eye and expose the lid margin. Rub the cleansing pad along the eyelid several times.3. Close the eye and rub pad vigorously over the eyelashes several times.4. Rinse thoroughly with clean, warm water and pat dry.5. Repeat the procedure for the other eye using a clean cleansing pad.
Availability:	<p>Lid-Care* is packaged in 100 mL bottle in a kit with 100 cleansing pads or boxes of 30 pre-moistened towelettes.</p>

Novartis Ophthalmics
Novartis Pharmaceuticals Canada Inc.
Mississauga, ON L5N 2X7

*Lid-Care is a registered trademark